

BACKGROUND

The formation of adhesions following abdominopelvic procedures is nearly unavoidable. Adhesions can lead to infertility and post-operative pain. Adept solution works as an instillate and forms a fluid reservoir in the peritoneal cavity with a prolonged residence time of up to 4 days and is slowly absorbed by the lymphatic system. Because of this, a large pivotal trial to confirm Adept's clinical efficacy and safety was performed.

STUDY OBJECTIVES

To evaluate the efficacy and safety of Adept in reducing adhesions after laparoscopic gynaecological surgery involving adhesiolysis while comparing it to Lactated Ringer's Solution (LRS).

STUDY DESIGN

Double-blind, randomised multicenter study consisting of four visits and conducted at 16 referral centers. Patients consisted of women ≥ 18 years old in general good health. Laparoscopic surgery was planned for a gynaecological procedure that included adhesiolysis followed by a second follow-up laparoscopy 4-8 weeks later. Details of each visit are as follows:

Visit 1 (Up to 4 weeks prior to surgery)	Visit 2 (First lap procedure)	Visit 3 (1-3 weeks after initial surgery)	Visit 4 (Final lap procedure 4-8 weeks after initial surgery)
Patient underwent physical exam and samples were taken for lab tests	Presence/absence of adhesions, extent & severity, and AFS* scores recorded	Patients underwent clinical laboratory tests and a physical examination	Videotaping and scoring of all available anatomical sites were performed, including AFS scores
Patient's medical history was recorded	Abdomen irrigated with minimum of 100 ml of study solution every 30 min.	Concomitant medications and any adverse events were recorded	All adverse events and concomitant medications were recorded
Degree of pain assessed at baseline using a visual analog scale (VAS)	A minimum of three adhesions had to be lysed and recorded at initial surgery		Degree of pain assessed at final visit using VAS
	Any remaining study solution was aspirated, and 1000 ml of study solution instilled		

Adhesions were scored at all 23 or all available anatomical sites. The presence or absence and severity of the adhesions were recorded. Extent was defined as one of the following:

Localized	Moderate	Extensive
< 1/3 of the adhesion site covered	1/3-2/3 of the adhesion site covered	>2/3 of the adhesion site covered

Severity was defined as:

Mild: filmy and avascular **Severe:** dense, cohesive, or vascular

Study efficacy measures included: clinical success, incidence, extent, and severity of adhesions and adhesion scoring using the AFS classification for adnexal adhesions. During a discussion with the FDA, the clinical success of adhesion reduction for a patient was determined as a reduction in adhesions of at least three or 30% of sites lysed (depending on which is greater) between initial surgery and the follow-up laparoscopy.

Safety was assessed by serious adverse events (SAEs), adverse events and changes in lab values. Patients also recorded their well-being and all concomitant meds.

RESULTS

The study began with 203 patients in the Adept group and 199 patients in the LRS group. More Adept patients achieved clinical success than did LRS patients (49% vs. 38%). Of the 57 patients treated with Adept and the 55 patients treated with LRS who showed moderate/severe AFS scores at initial surgery, 29 patients (51%) and 19 (35%) patients, respectively, moved to a minimal/mild AFS category due to improvements in their AFS score. Safety was comparable with both groups. The frequency of adverse events and the number of patients who reported them were similar in both treatment groups.

Patients treated with Adept had more favorable outcomes than those receiving LRS. In patients with infertility and in the subgroups with diagnoses of endometriosis or pain, the parameters of clinical success and AFS score also showed more favorable outcomes with Adept compared with LRS.


CLINICALLY PROVEN RESULTS

PAMELA study

Double-blind randomised controlled trial comparing ADEPT (n=203) and Lactated Ringer's Solution (n=199) during laparoscopic gynaecological

**ODDS RATIO
(p-value)**

12 out of 13 study parameters showed a favourability for ADEPT [Modified from Brown et al. 2007]



CONCLUSIONS

Based on this study, Adept is easy-to-apply, safe, and effective as an adhesion reduction agent in laparoscopy.

Using Adept as an irrigant and postoperative instillate reduces adhesions after laparoscopic gynaecologic adhesiolysis more than LRS. Although both Adept and LRS-treated patients showed a reduction in AFS score from their baseline score, the Adept patients displayed a more significant reduction in AFS score of 2.5-3.0 units beyond that of the LRS patients. A reduction of this amount is considered clinically significant since a patient can now be placed in a better prognostic category (minimal, mild) for pregnancy as opposed to what would have occurred otherwise (moderate, severe).

Any medical device product quality complaints (including medical device adverse incidents) relating to Baxter products can be reported directly to the Baxter Country Quality Assurance Team:

In the UK on +44 (0)1604 704603, or by email to UK_SHS_QA_Complaints@baxter.com

In Ireland on +353 (0)1 2065500 or by email to shs_complaints_dublin@baxter.com

Alternatively please report directly to your Baxter Representative, who will take the details and forward to the Baxter Country Quality Assurance Team.

Medical device adverse incidents should also be reported:

In the UK to the MHRA.

Reporting forms and information can be found at: www.mhra.gov.uk/safetyinformation/reportingsafetyproblems/index.htm

In Ireland to the HPRA.

Reporting forms and information can be found at: <http://www.hpra.ie/homepage/about-us/report-an-issue>