

ACTIFUSE Bone Graft Substitute Demonstrates Similar Fusion Rates to BMP-2 in Posterolateral Lumbar Fusion (PLF) Procedures

In a prospective, randomized, multicenter comparative study in 103 patients with degenerative spinal disorders requiring PLF (1-3 levels), ACTIFUSE Bone Graft Substitute provided fusion rates similar to BMP-2 (Level of Evidence:2).¹

Key Points

- No significant difference in successful fusion between ACTIFUSE Bone Graft Substitute and BMP-2 [INFUSE Bone Graft] at 12 and 24 months.
- Pain assessment (VAS Score) at 6 months was significantly reduced with ACTIFUSE vs. BMP-2.
- Successful neurological status was higher with ACTIFUSE vs. BMP-2 at 6 months.

Overview

Although autologous iliac crest bone has long been considered the gold standard bone graft in PLF, harvesting is associated with significant morbidity at the donor site and prolonged hospitalization. Hence, there is much interest in substitute biomaterials.^{2,3}

ACTIFUSE Bone Graft Substitute is a silicate substituted calcium phosphate synthetic bone graft that combines an optimized osteoconductive scaffold with an aqueous polymer gel resulting in a moldable, cohesive material that facilitates rapid and sustained bone ingrowth. Recombinant human bone morphogenetic protein, in combination with a bovine type I collagen scaffold, or BMP-2 [INFUSE, Medtronic, Inc.], is a bone inductive agent.

Vertebral fusion success rates for ACTIFUSE and BMP-2 were assessed by Coughlan et al¹ following PLF using computed tomography (CT) scans evaluated by independent radiographic reviewers at 12 and 24 months. CT scans were graded according to the method of Glassman et al.⁴ Clinical outcomes included pain scores on a visual analog scale (VAS), Oswestry Disability Index (ODI), SF-36 and adverse events.

Results

A total of 103 patients were enrolled and received treatment (ACTIFUSE, n=51; BMP-2, n=52); 96 patients completed the study. The per-protocol (PP) population, defined as all patients who completed the study without protocol deviations, comprised 62 patients with evaluable fusion results (ACTIFUSE, n=35; BMP-2, n=27) at the primary endpoint of 12-month follow-up. Baseline demographics were similar in both treatment groups.


The primary diagnoses were spondylolisthesis (n=49;47.6%) and degenerative disc disease (n=46;44.7%). For most patients (n=74;71.8%), there was one involved vertebral level; 28 (27.2%) patients had two levels and one patient (1.0%) had three levels (**Table 1**).

Table 1. Primary Diagnoses/Vertebral Levels

Diagnosis/Level	Number/Percentage
Spondylolisthesis	n=49/47.6%
Degenerative Disc Disease	n=46/44.7%
One Level	n=74/71.8%
Two Levels	n=28/27.2%
Three level	n=1/1%


At 12 months, in the PP population (62 patients), fusion success was achieved in 25 of 35 (71.4%) in the ACTIFUSE group and in 20 of 27 (74.1%) in the BMP-2 group (P=1.00, **Figure 1**).

Figure 1. Proportion of PP patients with fusion success at 12 and 24 months.



Representative CT images of solid fusion at 12 months are shown in **Figure 2**. At 24 months in the PP analysis, fusions were noted in 78.6% and 84.8% of patients in the ACTIFUSE and BMP-2 groups, respectively ($P=0.56$)(**Figure 1**).

Figure 2. CT scans at 12 months. (A) Patient treated with ACTIFUSE (i) central section (ii) outer section. (b) Patient treated with BMP-2 (i) central section (ii) outer section.



For both treatment groups, mean VAS scores for back pain were lower than preoperative VAS scores for all subsequent time-points (**Figure 3**); at 6 months, the ACTIFUSE VAS scores were significantly lower than the BMP-2 group. **Table 2** summarizes ODI scores in the two treatment groups throughout the study. In the ACTIFUSE and BMP-2 groups, ODI score and SF-36 showed steady improvements over time (SF-36 significant at 6 months). At all time-points, the proportion of patients deemed as experiencing neurological success was higher in the ACTIFUSE vs. BMP-2 group (**Figure 4**), although the difference was not significant, except for the 6-month timepoint. The numbers and frequencies of AEs were similar between the two treatment groups. The most frequently reported AE was pain, followed by wound secretion, pain in an extremity, back pain, nausea and procedural pain.

Figure 3. Mean visual analog scale (VAS) scores up to 24 months after treatment with ACTIFUSE or BMP-2.





Table 2. Oswestry Disability Index (ODI) Scores

VISIT	ACTIFUSE	BMP-2
Preoperative	48.5±16.6 (n=47)	48.5±16.1 (n=35)
6 weeks	41.8±15.9 (n=36)	47.1±17.6 (n=25)
3 months	33.7±18.0 (n=36)	36.3±20.9 (n=30)
6 months	28.3±17.2 (n=40)	36.1±20.5 (n=29)
12 months	29.4±18.8 (n=41)	36.6±22.6 (n=32)
24 months	29.1±21.0 (n=39)	29.6±21.0 (n=27)

Figure 4. Neurological status up to 24 months after treatment with ACTIFUSE or BMP-2.



Discussion and Conclusions

Although BMP-2 has not gained FDA approval for use in PLF, it has been demonstrated as highly efficacious at inducing bone ingrowth in humans.⁵⁻⁷ In patients undergoing PLF, equivalent bone fusion success was achieved using BMP-2 compared with control patients receiving iliac crest bone graft.⁸

The study was underpowered due to the number of patients who could not be included in the PP population. This was in part due to the number of protocol deviations, mostly arising from the omission of bulking agent to the BMP-2 material.

ACTIFUSE was safe and well tolerated in this study and provided fusion rates similar to those observed in patients receiving BMP-2 bone graft material. On the basis of historical control data, ACTIFUSE may be as useful as iliac crest bone autograft in the context of spine fusion surgery, with less risk of unwanted donor site morbidity.

References

1. Coughlan M, Davies M, Mostert AK, et al. A prospective, randomized, multicenter study comparing silicated calcium phosphate versus BMP-2 synthetic bone graft in posterolateral instrumented lumbar fusion for degenerative spine disorders. *Spine*. 2018;43:E860-868.
2. Miura Y, Imagama S, Yoda M, et al. Is local bone as viable as a source of bone graft in posterior lumbar interbody fusion? *Spine*. 2003;28:2386-2389.
3. Glassman SD, Howard JM, Sweet A, et al. Complications and concerns with osteobiologics for spine fusion in clinical practice. *Spine*. 2010;35:1621-1628.
4. Glassman SD, Dimar JR, Carreon LY, et al. Initial fusion rates with recombinant human bone morphogenetic protein-2/compression resistant matrix and a hydroxyapatite and tricalcium phosphate/collagen carrier in posterolateral spinal fusion. *Spine*. 2005;30:1694-1698
5. Boden SD, Kang J, Sandhu H, et al. Use of recombinant human bone morphogenetic protein-2 to achieve posterolateral lumbar spine fusion in humans: a prospective, randomized clinical pilot trial. 2002 Volvo Award in clinical studies. *Spine*. 2002;27:2662-2673.
6. Boden SD, Grob D, Damien C. Ne-Osteo bone growth factor for posterolateral lumbar spine fusion: results from a nonhuman primate study and a prospective human clinical pilot study. *Spine*. 2004;29:504-514.
7. Vaccaro AR, Patel T, Fischgrund J, et al. A pilot study evaluating safety and efficacy of OP-1 Putty (rhBMP-7) as a replacement for iliac crest autograft in posterolumbar arthrodesis for degenerative spondylolisthesis. *Spine*. 2004;29:1885-1892.
8. Glassman SD, Carrion L, Djurasovic M, et al. Posterolateral lumbar spine fusion with INFUSE bone graft. *Spine Journal*. 2007;7:44-49.

This study was funded by Baxter Healthcare.

Figures, tables and selected text were modified from Coughlan et al' under the Creative Commons Attribution 4.0 International License.

Baxter and Actifuse are registered trademarks of Baxter International Inc., or its subs
Infuse is a registered trademark of Medtronic Inc.

UKI/MG95/19-0006 December 2019

Any medical device product quality complaints (including medical device adverse incidents) relating to Baxter products can be reported directly to the Baxter Country Quality Assurance Team: In the UK on +44 (0)1604 704603, or by email to UK_SHS_QA_Complaints@baxter.com. In Ireland on +353 (0)1 2065500 or by email to shs_complaints_dublin@baxter.com

Alternatively please report directly to your Baxter Representative, who will take the details and forward to the Baxter Country Quality Assurance Team.

Medical device adverse incidents should also be reported: In the UK to the MHRA.

Reporting forms and information can be found at: www.mhra.gov.uk/safetyinformation/reportingsafetyproblems/index.htm . In Ireland to the HPRA. Reporting forms and information can be found at: <http://www.hpra.ie/homepage/about-us/report-an-issue>.

Baxter Healthcare Ltd
Unit 7 Deansgrange
Business park
Blackrock
Co. Dublin

Baxter Healthcare Ltd
Wallingford Road, Compton
Newbury, RG20 7QW
www.baxterhealthcare.co.uk

PharmedUK 
A Pharmed Group Company
8 Manor Park,
Wildmere Road,
Banbury, Oxfordshire, OX16 3TB,
United Kingdom

Baxter