

# Preparation Instructions (10 mL)

## FLOSEAL Haemostatic Matrix

# Floseal

## HAEMOSTATIC MATRIX

# 1

Open all packages and place the components in a **sterile field**.



Inspect the integrity of the contents of the FLOSEAL Kit. If the packaging or vial has been damaged or opened, do not use.



Open the package containing the Sodium Chloride ampoule and deliver the ampoule to the sterile field using an aseptic technique.



Open the outer pouch of the thrombin component and deliver the inner pouch to the sterile field. Once placed in the sterile field, the inner package may be opened at any time.



Open the outer pouch of the gelatin matrix component and deliver the inner pouch to the sterile field. Once placed in the sterile field, the inner package may be opened at any time.



Take care not to misplace the white luer connector.

Product	Ordering code
Floseal Haemostatic Matrix 10 mL	1503354
Laparoscopic Applicator, 6 pack	0600125

# 2

## Preparing the Thrombin Solution

Once mixed, the Thrombin Solution can be used immediately or may be stored at 2-25°C in the vial up to **4 hours**.

**Components:**



Remove the Tyvek lid from the vial adapter packaging.



Remove twist-off cap from the Sodium Chloride ampoule.



While gripping the vial adapter packaging, attach the vial adapter to the pre-filled Sodium Chloride ampoule using a luer lock motion and remove the vial adapter packaging.



Remove the plastic flip-off cap from the Thrombin vial.



While holding the vial adapter, pierce the rubber stopper of the Thrombin vial. Transfer the entire contents of the Sodium Chloride ampoule into the Thrombin vial.



Gently swirl the Thrombin vial with vial adapter and ampoule attached until the Thrombin is completely dissolved. Remove the empty Sodium Chloride ampoule and discard.



# 3

## Mixing the components

**Components:**



Attach the white luer connector to the empty 10mL syringe. With the white luer connector in place connect the empty syringe to the vial adapter on the Thrombin vial



Draw the thrombin solution from the vial to the indicated mark (8 mL). Discard the empty Thrombin vial and the vial adapter appropriately ensuring that the white luer connector remains connected to the syringe.



Remove the cap from the Gelatin Matrix Syringe carefully to avoid spilling the Gelatin Matrix granules. Connect the Gelatin Matrix Syringe to the syringe containing 8 mL of Thrombin solution using the white luer connector




Push the plunger of the Thrombin solution syringe to quickly pass the solution into the syringe containing the Gelatin Matrix. This constitutes "one pass." Transfer the Gelatin Matrix Thrombin solution mixture back and forth between the syringes for a total of 10 passes back and forth.



Ensure the syringe labeled FLOSEAL contains the reconstituted FLOSEAL Matrix. If desired, connect an Applicator Tip or Malleable Tip to the FLOSEAL syringe.



FLOSEAL may be used up to **8 hours** after mixing with the Thrombin Solution. Allow 30 seconds after preparation before FLOSEAL is applied to ensure optimal product consistency.



Any medical device product quality complaints (including medical device adverse incidents) relating to Baxter products can be reported directly to the Baxter Country Quality Assurance Team: In the UK on +44 (0)1604 704603, or by email to [UK\\_SHS\\_QA\\_Complaints@baxter.com](mailto:UK_SHS_QA_Complaints@baxter.com). In Ireland on +353 (0)1 2065500 or by email to [shs\\_complaints\\_dublin@baxter.com](mailto:shs_complaints_dublin@baxter.com). Alternatively please report directly to your Baxter Representative, who will take the details and forward to the Baxter Country Quality Assurance Team, Medical device adverse incidents should also be reported: In the UK to the MHRA. Reporting forms and information can be found at: [www.mhra.gov.uk/safetyinformation/reportingsafetyproblems/index.htm](http://www.mhra.gov.uk/safetyinformation/reportingsafetyproblems/index.htm). In Ireland to the HPRA. Reporting forms and information can be found at: <http://www.hpra.ie/homepage/about-us/report-an-issue>

UKI-AS19-230026 May 2023

Baxter Healthcare Ltd.  
Wallingford Road,  
Compton,  
Newbury, Berkshire,  
RG20 7QW

Baxter Healthcare Ltd.  
Unit 7,  
Deansgrange Business Park,  
Blackrock,  
Co Dublin

[www.baxterhealthcare.co.uk](http://www.baxterhealthcare.co.uk) [www.baxterhealthcare.ie](http://www.baxterhealthcare.ie)

**Baxter**