

Preparation Instructions (5 mL)

FLOSEAL Haemostatic Matrix

Floseal

HAEMOSTATIC MATRIX

1 Open all packages and place the components in a **sterile field**.



1 Inspect the integrity of the contents of the FLOSEAL Kit. If the packaging or vial has been damaged or opened, do not use.



2 Open the package containing the ampoule and deliver the ampoule to the sterile field.



3 Open the outer pouch of the thrombin component and deliver the inner pouch to the sterile field. Once placed in the sterile field, the inner package may be opened at any time.



4 Open the outer pouch of the gelatin matrix component and deliver the inner pouch to the sterile field. Once placed in the sterile field, the inner package may be opened at any time.

Product	Ordering code
Floseal Haemostatic Matrix 5 mL	1503353
Laparoscopic Applicator, 6 pack	0600125

2 Preparing the Thrombin Solution

Once mixed, the Thrombin Solution can be used immediately or may be stored at 2-25°C in the vial up to **4 hours**.

Components:



5 Remove the Tyvek lid from the vial adapter packaging.



6 Remove twist-off cap from the Sodium Chloride ampoule.



7 While gripping the vial adapter packaging, attach the pre-filled Sodium Chloride ampoule to the luer connector of the vial adapter and remove from remaining packaging.



8 Remove the plastic flip-off cap from the Thrombin vial.



9 While holding the vial adapter, pierce the rubber stopper of the Thrombin vial.

Transfer the entire contents of the Sodium Chloride ampoule into the Thrombin vial.



10 Gently swirl the Thrombin vial with vial adapter and ampoule attached until the Thrombin is completely dissolved.



Sodium Chloride ampoule

Thrombin vial

Needle-free vial adapter

3 Mixing the components

Components:



11 Remove Sodium Chloride ampoule.



12 Attach empty 5 mL syringe with female luer connector to the vial adapter and draw the thrombin solution from the vial to the indicated mark (4 mL).

Discard the Thrombin vial and the vial adapter appropriately.



13 Remove the luer cap from the FLOSEAL gelatin matrix syringe.

Connect this syringe to the syringe containing the Thrombin Solution.



14 Push the Thrombin Solution quickly into the gelatin matrix.

Transfer the mixture back and forth at least 10 times, making sure the mixed FLOSEAL material ends up in the FLOSEAL syringe.

Allow 30 seconds after preparation before FLOSEAL is applied to ensure optimal product consistency. FLOSEAL may be used up to **8 hours** after mixing with the Thrombin Solution.



15 FLOSEAL may be extruded directly from the syringe. If desired, connect an applicator tip to the FLOSEAL syringe.



FLOSEAL gelatin matrix syringe

Empty 5 mL syringe

Applicator tips

Malleable tip

Any medical device product quality complaints (including medical device adverse incidents) relating to Baxter products can be reported directly to the Baxter Country Quality Assurance Team: In the UK on +44 (0)1604 704603, or by email to UK_SHS_QA_Complaints@baxter.com. In Ireland on +353 (0)1 2065500 or by email to shs_complaints_dublin@baxter.com. Alternatively please report directly to your Baxter Representative, who will take the details and forward to the Baxter Country Quality Assurance Team.

Medical device adverse incidents should also be reported to:

In the UK to the MHRA, Reporting forms and information can be found at: www.mhra.gov.uk/safetyinformation/reportingsafetyproblems/index.htm.

In Ireland to the HPRA, Reporting forms and information can be found at: <http://www.hpra.ie/homepage/about-us/report-an-issue>

UKI-AS19-230025 May 2023

Baxter Healthcare Ltd.
Wallingford Road,
Compton,
Newbury, Berkshire,
RG20 7QW

www.baxterhealthcare.co.uk

Baxter Healthcare Ltd.
Unit 7,
Deansgrange Business Park,
Blackrock,
Co Dublin

www.baxterhealthcare.ie

Baxter