

Artiss

FIBRIN SEALANT (HUMAN)

Ready to use*

Preparation Guide

*3 Ways to Thaw the Frozen Pre-Filled PRIMA Syringe

*Plan Ahead: Thaw up to 14 days before use. Record the date and time on the sticker provided on the box.



Thawing / warming - ON sterile field:

Thawing/warming times 33°-37°C

- 2 mL in bath without pouch approximately 5 mins.
- 4 mL in bath without pouch approximately 5 mins.
- 10 mL in bath without pouch approximately 10 mins.

(Ensure the contents of the syringe remain completely submerged throughout thawing.)



Thawing / warming - OFF sterile field:

Thawing/warming times 33°-37°C

- 2 mL in bath with pouches approximately 15 mins.
- 4 mL in bath with pouches approximately 20 mins.
- 10 mL in bath with pouches approximately 35 mins.

(Ensure the pouches remain submerged throughout thawing.)



Thawing / warming - OFF sterile field (Incubator):


Thawing/warming times 33°-37°C

- 2 mL in pouches approximately 40 mins.
- 4 mL in pouches approximately 50 mins.
- 10 mL in pouches approximately 90 mins.

DO NOT microwave.


Application of ARTISS Fibrin Sealant [Human] must be completed within 4 hours of warming to 33°C-37°C. Discard any unused product.

Product Storage



Long term: Store at <-20°C

Short term: Unopened pouches, thawed at room temperature, may be stored for up to 14 days at room temperature (15°C to 25°C) after removal from the freezer. **Prior to application, the product must be warmed to 33-37°C.**



Once the package is opened or the product is warmed to 33°C to 37°C, **it must be used within 4 hours.**



DO NOT REFRIGERATE.

Once thawed, **DO NOT** refrigerate or re-freeze.

DO NOT expose to temperatures above 37°C.

DO NOT use after the expiration date.

